


# Is your hospital 100% safe?

## Identify all your MRSA-positive patients.

Abacus Diagnostica introduces a fully automated PCR-system for rapid MRSA and SA detection for hospitals that want to make an extra effort to stay free of MRSA.

We offer high sensitivity IVD testing on the fastest currently available instrument for screening and diagnosis.



### Various sample matrixes

- Multi and nasal swabs, blood culture, any culture plate
- 4 tests in 50 minutes


### Fast & Reliable

- Automatic result interpretation
- No need for expertise in PCR



### Safety = Peace of Mind

- Comprehensive strain coverage
- No false negatives



### Easy to use

- Preloaded dry-reagent PCR chemistry
- No risk of cross-contamination
- No need for dedicated premises

Assay steps for the GenomEra MRSA/SA Blood Culture test:



Mix the positive blood culture sample with buffer



Pipette to test chips



Start the GenomEra CDX assay run

Come and meet us in 2011 at

- IBMS, Birmingham, 26-28 Sep, stand # 840, Balcony level, hall 3.
- Near patient testing in microbiology, 10 Oct, London
- Infection Prevention Scottish National Conference, 27 Oct, Glasgow



Point Of Care Testing Ltd  
Phone: 01241 439 020  
Email: info@pocT.co.uk  
www.pocT.co.uk